



Sintec Optronics Pte Ltd

50 Bukit Batok Street 23 #07-15 Singapore 659578

Tel: +65 67781866 Fax: +65 65152908

Optical Tables

High-density plated steel honey comb structure is used and there is good vibration insulation. The optical tabletop consists of magnetic stainless steel and there is good flatness.

Specifications:

Tabletop thickness: 50mm

Flatness: 0.1mm/600mm

Mounting hole pattern: M6 on 25mm centers

Model	Tabletop size (mm)
OTBB33-1	300x300
OTBB36-1	300x600
OTBB39-1	300x900
OTBB312-1	300x1200
OTBB66-1	600x600
OTBB69-1	600x900
OTBB612-1	600x1200
OTBB99-1	900x900
OTBB912-1	900x1200
OTBB1212-1	1200x1200



OTBS69-1 Optical Table Frame

It is rigidly designed and provides a firm and stable work surface. The top of each leg is rubber covered to protect the bottom of the worktop. The height is adjustable using precise screws. It can support the optical tabletops with the sizes of 600x900mm, 900x900mm or 900x1200.

Specifications:

Normal height: 700mm

Height adjustment range: +30mm

-10mm

